


You have scheduled a visual rehabilitation session with the Eyesoft EMAA Home solution, please make sure you are familiar with the following steps before handing the headset to the patient.



WiFi Connection

To connect the headset to WiFi, choose "Configure", select your network, then enter the password using the central button on the right side of the headset. This step ensures that the Eyesoft solution is up to date.




Headset Login

Enter the patient code located on their file in your Eyesoft interface, then go to "Continue".



Scheduled Sessions

Once the patient's session is open in the headset, verify that their name is correct, then select "My sessions" to ensure that the scheduled sessions are present.



Time Limit

If the headset was not turned off correctly or remained on standby (see *On/Off button*), the following message may appear: *"Warning, you have exceeded the number of minutes... Please try again later."*

Turn off the headset then modify this limit from your interface, by accessing the "Configuration" tab on the patient file: adjust and save the number of minutes allowed per day.



Eyesoft support is available for orthoptists who have questions about EMAA Home. Requests sent directly by patients cannot be processed. Please read the user guide carefully before first use.